



### RNA-Seq Module 2: **FASTQ/SAM/BAM/BED file formats**

Arpad Danos, Felicia Gomez, Obi Griffith, Malachi Griffith, My Hoang, Mariam Khanfar, Chris Miller, Kartik Singhal Advanced Sequencing Technologies & Bioinformatics Analysis November 10-23, 2024








#### **Outline**


- What is a sequence read?
- Unaligned read
- Aligned read
- Difference between the stored aligned reads file formats

# Reads are the sequencer's best guess at what it saw for a given DNA molecule.


It's the "raw" data.



# Alignment is central to most genomics applications



# Alignment is central to most genomics applications



#### The FASTQ format

A "standard" format for storing and defining sequences from next-generation sequencing technologies.

```
Sequence ID
Sequence
Sequence
Sequence
Separator>
Quality scores
Sequence
SEQ_ID
GATTTGGGGTTCAAAGCAGTATCGATCAAATAGTAAATCCATTTGTTCAACTCACAGTTT

+
!''*((((***+))%%%++)(%%%%).1***-+*''))**55CCF>>>>>CCCCCCC65
```

- FASTQ files are generally used to store short-read data from high-throughput sequencing experiments
- The sequence and quality scores are usually put into a single line

### Sequence IDs


@HWUSI-EAS100R:6:73:941:1973#0/1

@EAS139:136:FC706VJ:2:2104:15343:197393 1:Y:18:ATCACG


| HWUSI-EAS100R | the unique instrument name  |
|---------------|---|
| 6             | flowcell lane   |
| 73            | tile number within the flowcell lane                                |
| 941           | 'x'-coordinate of the cluster within the tile                       |
| 1973          | 'y'-coordinate of the cluster within the tile                       |
| #0            | index number for a multiplexed sample (0 for no indexing)           |
| /1            | the member of a pair, /1 or /2 (paired-end or mate-pair reads only) |

| EAS139  | the unique instrument name   |
|---------|--|
| 136     | the run id   |
| FC706VJ | the flowcell id  |
| 2       | flowcell lane  |
| 2104    | tile number within the flowcell lane                                   |
| 15343   | 'x'-coordinate of the cluster within the tile                          |
| 197393  | 'y'-coordinate of the cluster within the tile                          |
| 1       | the member of a pair, 1 or 2 (paired-end or mate-pair reads only)      |
| Y       | Y if the read is filtered, N otherwise                                 |
| 18      | 0 when none of the control bits are on, otherwise it is an even number |
| ATCACG  | index sequence   |
|         |  |

### Quality scores




### Quality scores



ASCII is inferred by the probability of error ~ Q score (Phred score)

### Quality scores



Qualities are based on the Phred scale and are encoded

$$Q = -10*log_{10}(P_{err})$$
  
 $Q = -10*log_{10}(0.01)$   
 $Q = 20$  (Q is the Phred score)


- FASTQ files encodes phred scores as ASCII characters
- Phred quality scores characterize the quality of DNA sequences these scores are assigned by the sequencer
- A quality score of 20 (Q20) represents an error rate of 1 in 100 (meaning every 100 bp sequencing read may contain an error); call accuracy of 99% courtesy of Andrew Farrell

### Quality score encoding

| Dec | Hex | Char             | Dec        | Hex        | Char  | Dec  | Hex | Char | Dec | Hex | Char |
|-----|-----|------------------|------------|------------|-------|------|-----|------|-----|-----|------|
| 0   | 00  | Null             | 32         | 20         | Space | 64   | 40  | 0    | 96  | 60  | `    |
| 1   | 01  | Start of heading | I 33       | 21         | 1     | 65   | 41  | A    | 97  | 61  | a    |
| 2   | 02  | Start of text    | 34         | 22         |       | 66   | 42  | В    | 98  | 62  | b    |
| 3   | 03  | End of text      | 35         | 23         | #     | 67   | 43  | С    | 99  | 63  | c    |
| 4   | 04  | End of transmit  | 36         | 24         | ş     | 68   | 44  | D    | 100 | 64  | d    |
| 5   | 05  | Enquiry          | 37         | 25         | *     | I 69 | 45  | E    | 101 | 65  | e    |
| 6   | 06  | Acknowledge      | 38         | 26         | ٤     | 70   | 46  | F    | 102 | 66  | £    |
| 7   | 07  | Audible bell     | 39         | 27         | 1     | 71   | 47  | G    | 103 | 67  | g    |
| 8   | 08  | Backspace        | 40         | 28         | (     | 72   | 48  | H    | 104 | 68  | h    |
| 9   | 09  | Horizontal tab   | 41         | 29         | )     | I 73 | 49  | Ι    | 105 | 69  | i    |
| 10  | 0A  | Line feed        | 42         | 2A         | *     | 74   | 4A  | J    | 106 | 6A  | j    |
| 11  | OB  | Vertical tab     | <b>4</b> 3 | 2 B        | +     | 75   | 4B  | K    | 107 | 6B  | k    |
| 12  | OC. | Form feed        | 44         | 2 C        | ,     | 76   | 4C  | L    | 108 | 6C  | 1    |
| 13  | OD  | Carriage return  | 45         | 2 D        | -     | 77   | 4D  | M    | 109 | 6D  | m    |
| 14  | OE  | Shift out        | 46         | 2 <b>E</b> |       | 78   | 4E  | N    | 110 | 6E  | n    |
| 15  | OF  | Shift in         | 47         | 2 <b>F</b> | /     | 79   | 4F  | 0    | 111 | 6F  | 0    |
| 16  | 10  | Data link escape | 48         | 30         | 0     | 80   | 50  | P    | 112 | 70  | р    |
| 17  | 11  | Device control 1 | 49         | 31         | 1     | 81   | 51  | Q    | 113 | 71  | đ    |
| 18  | 12  | Device control 2 | 50         | 32         | 2     | 82   | 52  | R    | 114 | 72  | r    |
| 19  | 13  | Device control 3 | 51         | 33         | 3     | 83   | 53  | ន    | 115 | 73  | s    |
| 20  | 14  | Device control 4 | 52         | 34         | 4     | 84   | 54  | Т    | 116 | 74  | t    |
| 21  | 15  | Neg. acknowledge | 53         | 35         | 5     | 85   | 55  | U    | 117 | 75  | u    |
| 22  | 16  | Synchronous idle | 54         | 36         | 6     | 86   | 56  | V    | 118 | 76  | v    |
| 23  | 17  | End trans, block | 55         | 37         | 7     | 87   | 57  | W    | 119 | 77  | w    |
| 24  | 18  | Cancel           | 56         | 38         | 8     | 88   | 58  | X    | 120 | 78  | x    |
| 25  | 19  | End of medium    | 57         | 39         | 9     | 89   | 59  | Y    | 121 | 79  | У    |
| 26  | 1A  | Substitution     | 58         | 3A         | :     | 90   | 5A  | Z    | 122 | 7A  | z    |
| 27  | 1B  | Escape           | 59         | 3 B        | ;     | 91   | 5B  | [    | 123 | 7B  | {    |
| 28  | 1C  | File separator   | 60         | 3 C        | <     | 92   | 5C  | ١    | 124 | 7C  | I    |
| 29  | 1D  | Group separator  | 61         | 3 D        | =     | 93   | 5D  | ]    | 125 | 7D  | }    |
| 30  | 1E  | Record separator | 62         | 3 E        | >     | 94   | 5E  | ^    | 126 | 7E  | ~    |
| 31  | 1F  | Unit separator   | 63         | 3 <b>F</b> | ?     | 95   | 5F  | _    | 127 | 7F  |      |

Formula for getting PHRED quality from encoded quality:





10

- ASCII = American
  Standard Code for
  Information
  Interchange
- Every text symbol must have an integer value representing it inside the computer
- An ASCII code is the numerical representation of a character such as 'a' or '@'

Image courtesy of Andrew Farrell

### $Q=-10Log_{10}(P_{error})$

| Probability | of Error | Q  | <ul> <li>Higher Q scores indicate a<br/>smaller probability of</li> </ul> |
|-------------|----------|----|---|
| 1/1,000,000 | 0.00001  | 60 | error.  |
| 1/100,000   | 0.000010 | 50 | • Lower Q scores indicate   |
| 1/10,000    | 0.000100 | 40 | lower confidence in the called base.                                      |
| 1/1,000     | 0.001000 | 30 | <ul> <li>Increased false-positive variant calls</li> </ul>                |
| 1/100       | 0.010000 | 20 | <ul> <li>Q30 is a standard a<br/>benchmark for quality in</li> </ul>      |
| 1/10        | 0.100000 | 10 | benchmark for quality in next-generation sequencing                       |
| 1/1         | 1.000000 | 0  |   |

#### **FASTQ Report Summary**

**№**FastQC Report

Tue 12 Jan 2016 2009-08\_lib324\_miseq\_r0030\_251bp\_R1.fastq.gz

#### **Summary**

Basic Statistics

Per base sequence quality

Per tile sequence quality

Per sequence quality scores

Per base sequence content

Per sequence GC content

Per base N content

Sequence Length Distribution

Sequence Duplication Levels

Overrepresented sequences


Adapter Content

Kmer Content

#### Basic Statistics

| Measure                           | Value  |
|-----------------------------------|--|
| Filename                          | 2009-08_lib324_miseq_r0030_251bp_R1.fastq.gz |
| File type                         | Conventional base calls                      |
| Encoding                          | Sanger / Illumina 1.9                        |
| Total Sequences                   | 503810                                       |
| Sequences flagged as poor quality | 0  |
| Sequence length                   | 35-251                                       |
| %GC                               | 66   |

#### ②Per base sequence quality




#### FASTA format

We start with a reference genome to map to

The reference sequence (chromosome)

Sequence description



http://en.wikipedia.org/wiki/FASTA\_format


# Aligning to a reference genome; the crucial first step





Could fit here - but there are differences


Could fit here as well.


### Single-end alignment


















One mate maps uniquely, the other is unmapped



#### What needs to be stored?



#### Where did the read map? How confident are we that we are correct?

Which strand does the read come from?

Are there any differences with the reference?

What is the DNA sequence?

What are the quality scores for each base in the read?

What do we know about the mate?

Which read group does the read belong to?

http://samtools.github.io/hts-specs/SAMv1.pdf

#### What needs to be stored?



Where did the read map?
How confident are we that we are correct?

#### Which strand does the read come from?

Are there any differences with the reference?

What is the DNA sequence?

What are the quality scores for each base in the read?

What do we know about the mate?

Which read group does the read belong to?

http://samtools.github.io/hts-specs/SAMv1.pdf

Slide courtesy of Andrew Farrell

25

#### What needs to be stored?


Where did the read map?
How confident are we that we are correct? Which strand does the read come from?

Are there any differences with the reference?
What is the DNA sequence?
What are the quality scores for each base in the read?
What do we know about the mate?
Which read group does the read belong to?

26

http://samtools.github.io/hts-specs/SAMv1.pdf

### Store the alignment



#### **Standardize alignment formats**

SAM : Sequence Alignment/Map

BAM: Binary Alignment/Map

- Can be indexed allowing fast access of regions
- Simple format

27

- Can represent single and paired end reads
- Many toolkits now available to process data

http://samtools.github.io/hts-specs/SAMv1.pdf

#### **SAM/BAM format**

| Aspect       | SAM                         | BAM                              |
|--------------|-----------------------------|----------------------------------|
| Format       | Plain text (human-readable) | Binary (compressed)              |
| File Size    | Larger, less efficient      | Smaller, more efficient          |
| Access Speed | Slower for large files      | Faster, especially with indexing |
| Readability  | Directly readable           | Needs tools (e.g., samtools )    |

- BAM is a compressed version of SAM; lossless BGZF format
- BAM files are usually 'indexed'
  - A '.bai' file will be found beside the '.bam' file
  - Indexing provides fast retrieval of alignments overlapping a specified region without going through all alignments.
  - BAM must be sorted by the reference ID and then the leftmost coordinate before indexing

http://samtools.sourceforge.net/SAM1.pdf

#### **Example of SAM/BAM file format**

#### Example SAM/BAM/CRAM header section (abbreviated)

```
mgriffit@linus270 >> samtools view -H /qscmnt/qc13001/info/model data/2891632684/build136494552/alignments/136080019.bam | grep -P "SN\:22|HD|RG|PG"
        VN:1.4 SO:coordinate
                                UR:ftp://ftp.ncbi.nih.gov/genbank/genomes/Eukaryotes/vertebrates_mammals/Homo_sapiens/GRCh37/special_requests/GRCh37-lite.fa.gz AS:GRCh37-lite M5:a718acaa6135fdca8357d5bfe9
        SN:22 LN:51304566
4211dd SP:Homo sapiens
                                                                                                   PI:365 DS:paired end DT:2012-10-03T19:00:00-0500
       ID:2888721359
                        PL:illumina
                                         PU:D1BA4ACXX.3 LB:H KA-452198-0817007-cDNA-3-lib1
                                                                                                                                                             SM:H_KA-452198-0817007 CN:WUGSC
        ID: 2888721359
                       VN:2.0.8
                                         CL:tophat --library-type fr-secondstrand --bowtie-version=2.1.0
       ID:MarkDuplicates
                                                         PP:2888721359 WN:1.85(exported)
                                 PN:MarkDuplicates
                                                                                                   CL:net.sf.picard.sam.MarkDuplicates INPUT=[/gscmnt/gc13001/info/build merged alignments/merged-alignment-blad
e10-2-5.gsc.wustl.edu-jwalker-15434-136080019/scratch-ILg6Y/H_KA-452198-0817007-cDNA-3-lib1-2888360300.bam] OUTPUT=/gscmnt/gc13001/info/build_merged_alignments/merged-alignment-blade10-2-5.gsc.wustl.edu-jwalker-15434-136080019/scratch-ILg6Y/H_KA-452198-0817007-cDNA-3-lib1-2888360300.bam]
alker-15434-136080019/scratch-ILq6Y/H_KA-452198-0817007-cDNA-3-lib1-2888360300-post_dup.bam METRICS_FILE=/gscmnt/gc13001/info/build_merged_alignments/merged-alignment-blade10-2-5.gsc.wustl.edu-jwalker-1543
4-136080019/staging-liuJS/H KA-452198-0817007-cDNA-3-lib1-2888360300.metrics REMOVE DUPLICATES=false ASSUME SORTED=true MAX FILE HANDLES FOR READ ENDS MAP=9500 TMP DIR=[/gscmnt/gc13001/info/build merged al
ignments/merged-alignment-blade10-2-5.gsc.wustl.edu-jwalker-15434-136080019/scratch-ILg6Y] VALIDATION STRINGENCY=SILENT MAX RECORDS IN RAM=500000 PROGRAM RECORD ID=MarkDuplicates PROGRAM GROUP NAME=Mark
Duplicates MAX SEQUENCES FOR DISK READ ENDS MAP=50000 SORTING COLLECTION SIZE RATIO=0.25 READ NAME REGEX=[a-zA-Z0-9]+:[0-9]+):([0-9]+):([0-9]+):([0-9]+):* OPTICAL DUPLICATE PIXEL DISTANCE=100 VERBOSITY=INFO
QUIET=false COMPRESSION LEVEL=5 CREATE INDEX=false CREATE MD5 FILE=false
mgriffit@linus270 ~>
```

#### Example SAM/BAM/CRAM alignment section (only 10 alignments shown)

```
mgriffit@linus270 ->> samtools view -f 3 -F 1804 /gscmnt/gc13001/info/model_data/2891632684/build136494552/alignments/136080019.bam | head
HWI-ST495 129147882:3:2114:15769:38646 99
                                           11306 3
                                                       100M
                                                                               ACTGCGGGGCCCTCTTGCTTACTGTATAGTGGTGGCACGCCGCCTGCTGGCAGCTAGGGACATTGCAGGGTCCTCTTGCTCAAGGTGTAGTGGCAGCACGC
CC: Z: 15 MD: Z: 5A94
                                                                                                  PG:Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
      XN:i:0 X0:i:0 CP:i:102519765 AS:i:-5 XS:A:+ YT:Z:UU
HWI-ST495_129147882:3:2114:15769:38646 147
                                           11508 3
                                                                   11306
                                                                               ;5:CDCDCDECEFCD@9E=?7EEIIIIHCEGGIJJJJIIJJIHF@?00IHHFFGG?*JJJIJGHGEIJJIJJJJJJIHHCIEJJJHFHHGHFFEDFCCB
                                                                               CC:Z:15 MD:Z:34A65
                                                                                                  PG: Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
      XN:i:0 X0:i:0 CP:i:102519563 AS:i:-6 XS:A:+ YT:Z:UU
HWI-ST495 129147882:3:1210:1257:16203 163
                                          11810 3
                                                                               CCFFFFFHFHAFGGIIIJJJEEHGIGGGIJIJJGI?@EHIGIJDGHIHIGGIJJJJJJJJJJJJJGHHHGHFFFCDDDDDDCDCCCCCA;>@>@AA@:AA>AA
                                                                               CC: Z: 15 MD: Z: 100
                                                                                                  PG:Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:102519261 AS:i:0 XS:A:- YT:Z:UU
HWI-ST495 129147882:3:1210:1257:16203 83
                                          12055 3
                                                                               GAGCACTGGAGTGGAGTTTTCCTGTGGAGAGGAGCCATGCCTAGAGTGGGATGGGCCCATTGTTCATCTTCTGGCCCCCTGTTGTCTGCATGTAACTTAATAC
                                   1
                                                                   11810
CC:Z:15 MD:Z:100
                                                                                                  PG:Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:102519016 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495_129147882:3:2111:3117:78828 163
                                           12634 3
                                                       100M
                                                                   12746
                                                                               CC:Z:15 MD:Z:85G14
                                                                                                  PG:Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
      XN:i:0 X0:i:0 CP:i:102518437 AS:i:-5 XS:A:- YT:Z:UU
                                                                               {\tt GGGAGTGGCGTCGCCCCTAGGGCTCTACGGGGCCGGCATCTCCTGTCTCCTGGAGAGGCTTCGATGCCCCTCCACACCCTCTTGATCTTCCCTGTGATGTD}
HWI-ST495_129147882:3:2111:3117:78828 83
                                           12746 3
                                                                   12634 -212
DCABDBDDDDDDDDDDDDDDDDBDB@BDDDB@; CCCCCDEFD@; .?<HIGGEIGEHIGJJJJIIGIGIIHEGFEHFJIIIIIGJJJJHHHHHFFFFFC@
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
                                                                               CC:Z:15 MD:Z:37G62
                                                                                                  PG:Z:MarkDuplicates
      XN:i:0 X0:i:0 CP:i:102518325 AS:i:-5 XS:A:- YT:Z:UU
HWI-ST495_129147882:3:1102:4242:26638 99
                                          13503 3
                                                                   13779 376
                                                                               CGCTGTGCCCTTCCTTTGCTCTGCCCGCTGGAGACGGTGTTTGTCATGGGCCTGGTCTGCAGGGATCCTGCTACAAAGGTGAAACCCAGGAGAGTGTGGAC
CC: Z:2 MD: Z: 100
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
                                                                                                  PG:Z:MarkDuplicates
     XN:i:0 X0:i:0 CP:i:114357414 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495_129147882:3:1309:15328:74082 99
                                          13534 3
                                                                   13780
                                                                               AGACGGTGTTTGTCATGGGCCTGGTCTGCAGGGATCCTGCTACAAAGGTGAAACCCAGGAGAGTGTGGAGTCCAGAGTGTTGCCAGGACCCAGGCACAGG@
CC: Z: 2 MD: Z: 100
                                                                                                  PG:Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:114357383 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495_129147882:3:1308:10126:19636 99
                                           13779 3
                                                                   14027
                                                                               CCTCTGCAGGAGGCTGCCATTTGTCCTGCCCACCTTCTTAGAAGCGAGACGGAGCCCATCTGCTACTGCCCTTTCTATAATAACTAAAGTTAGCTGC
CC:Z:2 MD:Z:100
                                                                                                  PG:Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:114357140 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495 129147882:3:1102:4242:26638 147
                                                                               CCTCTGCAGGAGGCTGCCATTTGTCCTGCCCACCTTCTTAGAAGCGAGACGGAGCAGACCCATCTGCTACTGCCCCTTTCTATAATAACTAAAGTTAGCTG#
                                          13779 3
                                                                   13503
##DCCDDDCCBBBABCCDDDCBDDBBDHC?=GIIJIIIIJJGIIIIJJHJJJJJJGCIIJJJJJJIGHGJJIJJJJJJJJIIIIGGFGHHHHFFFFFCCC
                                                                               CC: Z:2 MD: Z: 100
                                                                                                  PG:Z:MarkDuplicates
                                                                                                                    RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:114357140 AS:i:0 XS:A:+ YT:Z:UU
mgriffit@linus270 ~>
```

## A SAM/BAM file is divided in header & alignment sections Example SAM/BAM header section (abbreviated)

```
mgriffit@linus270 >> samtools view -H /gscmnt/gc13001/info/model data/2891632684/build136494552/alignments/136080019.bam | grep -P "SN\:22|HD|RG|PG"
       VN:1.4 SO:coordinate
       SN:22 LN:51304566
                             UR:ftp://ftp.ncbi.nih.gov/genbank/genomes/Eukaryotes/vertebrates_mammals/Homo_sapiens/GRCh37/special_requests/GRCh37-lite.fa.gz AS:GRCh37-lite M5:a718acaa6135fdca8357d5bfe9
4211dd SP: Homo sapiens
       ID:2888721359 PL:illumina
                                                                                         PI:365 DS:paired end
                                     PU:D1BA4ACXX.3 LB:H_KA-452198-0817007-cDNA-3-lib1
                                                                                                               DT:2012-10-03T19:00:00-0500
                                                                                                                                             SM:H KA-452198-0817007 CN:WUGSC
      ID:2888721359 VN:2.0.8
                                     CL:tophat --library-type fr-secondstrand --bowtie-version=2.1.0
                                                                                         CL:net.sf.picard.sam.MarkDuplicates INPUT=[/gscmnt/gc13001/info/build_merged_alignments/merged-alignment-blad
       ID:MarkDuplicates
                             PN:MarkDuplicates
                                                    PP:2888721359 WN:1.85(exported)
e10-2-5.gsc.wustl.edu-jwalker-15434-136080019/scratch-ILg6Y/H_KA-452198-0817007-cDNA-3-lib1-2888360300.bam] OUTPUT=/gscmnt/gc13001/info/build_merged_alignments/merged-alignment-blade10-2-5.gsc.wustl.edu-jw
alker-15434-136080019/scratch-ILg6Y/H_KA-452198-0817007-cDNA-3-lib1-2888360300-post_dup.bam METRICS_FILE=/gscmnt/gc13001/info/build_merged_alignments/merged-alignment-blade10-2-5.gsc.wustl.edu-jwalker-1543
4-136080019/staging-liuJS/H_KA-452198-0817007-cDNA-3-lib1-2888360300.metrics REMOVE_DUPLICATES=false ASSUME_SORTED=true MAX_FILE_HANDLES_FOR_READ_ENDS_MAP=9500 TMP_DIR=[/gscmnt/gc13001/info/build_merged_al
ignments/merged-alignment-blade10-2-5.gsc.wustl.edu-jwalker-15434-136080019/scratch-ILg6Y] VALIDATION STRINGENCY=SILENT MAX RECORDS IN RAM=5000000 PROGRAM RECORD ID=MarkDuplicates PROGRAM GROUP NAME=Mark
Duplicates MAX_SEQUENCES_FOR_DISK_READ_ENDS_MAP=50000 SORTING_COLLECTION_SIZE_RATIO=0.25 READ_NAME_REGEX=[a-zA-Z0-9]+:[0-9]+):([0-9]+):([0-9]+):([0-9]+).* OPTICAL_DUPLICATE_PIXEL_DISTANCE=100 VERBOSITY=INFO
QUIET=false COMPRESSION LEVEL=5 CREATE INDEX=false CREATE MD5 FILE=false
mgriffit@linus270 ~>
                                          Version (VN) and sort order
                                          (SO) - Important!
                                                                                                                                   Reference sequence (SQ)
                                                                                                                                   and sequence length (LN)
                      @HD
                                              SO:coordinate
                      @SQ
                                  SN:20
                                             LN:63025520
                      @RG
                                  ID: HG00096
                                                         SM: HG00096
                                  ID: HG00096
                                                                     CL:/Users/AlistairNWard/Work/gkno/gkno launcher/tools/bwa/bwa mem -t
                          Read group (RG) and sample
                                                                      Programs (PG) that have
                           (SM)
                                                                      been run on the data
                                                                                                                                                               Slide courtesy of Andrew Farrell
```

#### SAM/BAM header section

- Used to describe source of data, reference sequence, method of alignment, etc.
- Each section begins with character '@' followed by a two-letter record type code. These are followed by two-letter tags and values:
  - @HD The header line
    - VN: format version
    - SO: Sorting order of alignments
  - @SQ Reference sequence dictionary
    - SN: reference sequence name
    - LN: reference sequence length
    - SP: species

- @RG Read group
  - ID: read group identifier
  - CN: name of sequencing center
  - SM: sample name
- @PG Program
  - PN: program name
  - VN: program version

#### A BAM file is divided in header and alignment sections Example SAM/BAM alignment section (only 10 alignments shown)

```
mgriffit@linus270 >> samtools view -f 3 -F 1804 /gscmnt/gc13001/info/model_data/2891632684/build136494552/alignments/136080019.bam | head
HWI-ST495 129147882:3:2114:15769:38646 99
                                             11306 3
                                                          100M =
                                                                       11508 302
                                                                                    ACTGCGGGGCCCTCTTGCTTACTGTATAGTGGTGGCACGCCGCCTGCTGGCAGCTAGGGACATTGCAGGGTCCTCTTGCTCAAGGTGTAGTGGCAGCACGC
CC: Z: 15 MD: Z: 5A94
                                                                                                       PG: Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
      XN:i:0 X0:i:0 CP:i:102519765 AS:i:-5 XS:A:+ YT:Z:UU
HWI-ST495_129147882:3:2114:15769:38646 147
                                             11508 3
                                                                                    -302
;5:CDCDCDECEFCD@9E=?7EEIIIIHCEGGIJJJJIIJJIHF@?00IHHFFGG?*JJJIJJGHGEIJJIJJJJJJIHHCIEJJJHFHHGHFFEDFCCB
                                                                                    CC:Z:15 MD:Z:34A65
                                                                                                       PG:Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
      XN:i:0 X0:i:0 CP:i:102519563 AS:i:-6 XS:A:+ YT:Z:UU
HWI-ST495 129147882:3:1210:1257:16203 163
                                             11810 3
                                                          100M
                                                                       12055
                                                                                    345
CCFFFFFHFHAFGGIIIJJJEEHGIGGGIJIJJGI?@EHIGIJDGHIHIGGIJJJJJJJJJJJJJGHHHGHFFFCDDDDDDCDCCCCCA;>@>@AA@:AA>AA
                                                                                    CC: Z: 15 MD: Z: 100
                                                                                                       PG:Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:102519261 AS:i:0 XS:A:- YT:Z:UU
HWI-ST495 129147882:3:1210:1257:16203 83
                                             12055 3
                                                                                    GAGCACTGGAGTGGAGTTTTCCTGTGGAGAGGAGCCATGCCTAGAGTGGGATGGGCCATTGTTCATCTTCTGGCCCCCTGTTGTCTGCATGTAACTTAATAC
CC: Z: 15 MD: Z: 100
                                                                                                       PG:Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:102519016 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495 129147882:3:2111:3117:78828
                                             12634 3
                                                                            212
                                                                                    163
@GFFFFFDHHHH9FHGIIFGAFDHEGII>GHIIIIIIIIIIIIIIIIIIIFHDDFFEEECEECCCACCCCC:AADCCBCC>CAC<CCCCC:@CB@@BAB##
                                                                                    CC: Z:15 MD: Z:85G14
                                                                                                       PG:Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
      XN:i:0 X0:i:0 CP:i:102518437 AS:i:-5 XS:A:- YT:Z:UU
HWI-ST495 129147882:3:2111:3117:78828 83
                                             12746 3
                                                                       12634
                                                                             -212
                                                                                    GGGAGTGGCGTCGCCCTAGGGCTCTACGGGGCCGGCATCTCCTGTCTCCTGGAGAGGCTTCGATGCCCCTCCACACCCTCTTGATCTTCCCTGTGATGTD
DCABDBDDDDDDDDDDDDDDDDBDB@BDDDB@;CCCCCDEFD@;.?<HIGGEIGEHIGJJJIIGIGIIHEGFEHFJIIIIIGJJJJHHHHHFFFFFC@@
                                                                                    CC: Z:15 MD: Z:37G62
                                                                                                       PG:Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:1 XM:i:
      XN:i:0 X0:i:0 CP:i:102518325 AS:i:-5 XS:A:- YT:Z:UU
HWI-ST495 129147882:3:1102:4242:26638 99
                                             13503 3
                                                                                    CGCTGTGCCCTTCCTTTGCTCTGCCCGCTGGAGACGGTGTTTGTCATGGGCCTGGTCTGCAGGGGATCCTGCTACAAAGGTGAAACCCAGGAGAGTGTGGAC
CCFFFFFHHHHHJJJIJJJJJJJJJJJJJJJGIIIIJJFHGGIJGIJJJJEGIJIJJHHIHHGHFFEFDEEECCCAACDDACDCDDDDDB?8?<B>A@CDC
                                                                                    CC: Z:2 MD: Z:100
                                                                                                       PG: Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:114357414 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495 129147882:3:1309:15328:74082 99
                                             13534 3
                                                                                    AGACGGTGTTTGTCATGGGCCTGGTCTGCAGGGATCCTGCTACAAAGGTGAAACCCAGGAGAGTGTGGAGTCCAGAGTGTTGCCAGGACCCAGGCACAGG
CC:Z:2 MD:Z:100
                                                                                                       PG: Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:114357383 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495 129147882:3:1308:10126:19636 99
                                                                             348
                                                                                    CCTCTGCAGGAGGCTGCCATTTGTCCTGCCCACCTTCTTAGAAGCGAGACGGAGCAGACCCATCTGCTACTGCCCTTTCTATAATAACTAAAGTTAGCTGC
CC:Z:2 MD:Z:100
                                                                                                       PG: Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:114357140 AS:i:0 XS:A:+ YT:Z:UU
HWI-ST495 129147882:3:1102:4242:26638
                               147
                                             13779 3
                                                                                    CCTCTGCAGGAGGCTGCCATTTGTCCTGCCCACCTTCTTAGAAGCGAGACGGAGCAGACCCATCTGCTACTGCCCTTTCTATAATAACTAAAGTTAGCTG#
##DCCDDDCCBBBABCCDDDCBDDBBDHC?=GIIJIIIIJIGIIIIJJHJJIJJIGCIIJJJJJJIGHGJJIJJJJJJJIIIIGGFGHHHHFFFFFCCC
                                                                                    CC:Z:2 MD:Z:100
                                                                                                       PG:Z:MarkDuplicates
                                                                                                                           RG:Z:2888721359 XG:i:0 NH:i:2 HI:i:0 NM:i:0 XM:i:
      XN:i:0 X0:i:0 CP:i:114357140 AS:i:0 XS:A:+ YT:Z:UU
mgriffit@linus270 ~>
                                                                                                                            Slide courtesy of Obi and Malachi Griffith
```

### SAM/BAM alignment section

| Col | Field        | Type                 | Regexp/Range            | Brief description                     |
|-----|--------------|----------------------|-------------------------|---------------------------------------|
| 1   | QNAME        | String               | [!-?A-~]{1,255}         | Query template NAME                   |
| 2   | FLAG         | $\operatorname{Int}$ | $[0,2^{16}-1]$          | bitwise FLAG                          |
| 3   | RNAME        | String               | \* [!-()+-<>-~][!-~]*   | Reference sequence NAME               |
| 4   | POS          | Int                  | $[0,2^{29}-1]$          | 1-based leftmost mapping POSition     |
| 5   | MAPQ         | Int                  | $[0,2^8-1]$             | MAPping Quality                       |
| 6   | CIGAR        | String               | \* ([0-9]+[MIDNSHPX=])+ | CIGAR string                          |
| 7   | <b>RNEXT</b> | String               | \* = [!-()+-<>-~][!-~]* | Ref. name of the mate/next segment    |
| 8   | <b>PNEXT</b> | $\mathbf{Int}$       | $[0,2^{29}-1]$          | Position of the mate/next segment     |
| 9   | TLEN         | $\operatorname{Int}$ | $[-2^{29}+1,2^{29}-1]$  | observed Template LENgth              |
| 10  | SEQ          | String               | \* [A-Za-z=.]+          | segment SEQuence                      |
| 11  | QUAL         | String               | [!-~]+                  | ASCII of Phred-scaled base QUALity+33 |

| RR062634 | .14576120 | 163           | 20       | 899919               | 60     | 100M              | =     | 900037       | 218     | TTCCCCAGTAGCTGGGATTACAGGCATACGCCACC  |
|----------|-----------|---------------|----------|----------------------|--------|-------------------|-------|--------------|---------|--------------------------------------|
| 1        | QNAME     | e.g. HWI-ST49 | 95 1291  | 47882 <b>:</b> 1:230 | 2:102  | 69 <b>:</b> 12362 |       |              |         |                                      |
| 2        | FLAG      | e.g. 163      | _        |                      |        |                   |       |              |         |                                      |
| 3        | RNAME     | e.g. 20       |          |                      |        |                   |       |              |         |                                      |
| 4        | POS       | e.g. 899919   |          |                      |        |                   |       |              |         |                                      |
| 5        | MAPQ      | e.g. 60       |          |                      |        |                   |       |              |         |                                      |
| 6        | CIGAR     | e.g. 100M     |          |                      |        |                   |       |              |         |                                      |
| 7        | RNEXT     | e.g. =        |          |                      |        |                   |       |              |         |                                      |
| 8        | PNEXT     | e.g. 900037   |          |                      |        |                   |       |              |         |                                      |
| 9        | TLEN      | e.g. 218      |          |                      |        |                   |       |              |         |                                      |
| 10       | SEQ       | e.g. CCTGTTTC | CTCCACA  | AAGTGTTTACT          | TTTGG  | ATTTTTGCCAG       | TCTA  | CAGGTGAAGCCC | TGGAGA' | TTCTTATTAGTGATTTGGGCTGGGGCCTGGCCATGT |
| 11       | QUAL      | e.g. CCCFFFF  | ГНННННЈ. | JIJFIJJJJJJ          | JJJJJI | HIJJJJJJJIJ       | JJJJG | GHIJHIJJJJJJ | JJJGHG  | GIJJJJJJIJEEHHHHFFFFCDCDDDDDDDB@ACDD |

#### SAM/BAM flags explained

- 12 bitwise flags describing the alignment
- Stored as a binary string of length 12 instead of 12 columns of data
- Value of '1' indicates the flag is set. e.g. 001000000000
- All combinations can be represented as a number from 0 to 4095 (i.e. 2<sup>12</sup>-1). This number is used in the BAM/SAM file.
- You can specify 'required' or 'filter' flags in samtools view using the '-f' and '-F' options respectively

| Bit  |       | Description  |  |  |  |  |  |
|------|-------|--|--|--|--|--|--|
| 1    | 0x1   | template having multiple segments in sequencing                    |  |  |  |  |  |
| 2    | 0x2   | each segment properly aligned according to the aligner             |  |  |  |  |  |
| 4    | 0x4   | segment unmapped   |  |  |  |  |  |
| 8    | 0x8   | next segment in the template unmapped                              |  |  |  |  |  |
| 16   | 0x10  | SEQ being reverse complemented                                     |  |  |  |  |  |
| 32   | 0x20  | SEQ of the next segment in the template being reverse complemented |  |  |  |  |  |
| 64   | 0x40  | the first segment in the template                                  |  |  |  |  |  |
| 128  | 0x80  | the last segment in the template                                   |  |  |  |  |  |
| 256  | 0x100 | secondary alignment  |  |  |  |  |  |
| 512  | 0x200 | not passing filters, such as platform/vendor quality controls      |  |  |  |  |  |
| 1024 | 0x400 | PCR or optical duplicate   |  |  |  |  |  |
| 2048 | 0x800 | supplementary alignment  |  |  |  |  |  |

Note that to maximize confusion, each bit is described in the SAM specification using its hexadecimal representation (i.e., 0x10' = 16 and 0x40' = 64).

http://broadinstitute.github.io/picard/explain-flags.html

#### **CIGAR** strings explained

- •The 'CIGAR' (**C**ompact **I**diosyncratic **G**apped **A**lignment **R**eport)
- •The CIGAR string is a sequence of base lengths and associated 'operations' indicating which bases align to the reference (either a match or mismatch), are deleted, are inserted, represent introns, etc.

| Op   | BAM      | Description   |  |  |  |  |  |  |  |
|------|----------|---|--|--|--|--|--|--|--|
| M    | 0        | alignment match (can be a sequence match or mismatch) |  |  |  |  |  |  |  |
| I    | 1        | insertion to the reference                            |  |  |  |  |  |  |  |
| D    | <b>2</b> | deletion from the reference                           |  |  |  |  |  |  |  |
| N    | 3        | skipped region from the reference                     |  |  |  |  |  |  |  |
| S    | 4        | soft clipping (clipped sequences present in SEQ)      |  |  |  |  |  |  |  |
| H    | 5        | hard clipping (clipped sequences NOT present in SEQ)  |  |  |  |  |  |  |  |
| P    | 6        | padding (silent deletion from padded reference)       |  |  |  |  |  |  |  |
| =    | 7        | sequence match  |  |  |  |  |  |  |  |
| X    | 8        | sequence mismatch                                     |  |  |  |  |  |  |  |
| Refe | erence:  | ACTTTTCATCCCTAAACAACCCTGTGTTTTCCC                     |  |  |  |  |  |  |  |
| San  | nple:    | ACCGGTCATTAAATTTCAACCCTTCTGTGAAATCCC                  |  |  |  |  |  |  |  |
|      |          | 9M3D4M3I5M3I12  |  |  |  |  |  |  |  |

rence: ACTTTTCATCCCTAAA---CAACC---CTGTGTTTTCCC
ple: ACCGGTCAT---TAAATTTCAACCCTTCTGTGAAATCCC
9M3D4M3I5M3I12M
2M3X4M

36

#### **CRAM** files

- CRAM is an ultra-compressed version of a BAM

   Usually between 30-60% smaller than the
   corresponding BAM
- Stores "diffs" from the reference genome
  - requires the matching reference genome to restore original data!
- Base quality binning may be used as well
- Some tools still require conversion back to bam

#### SAM > BAM > CRAM

| Quality Score Bins | Example of Empirically<br>Mapped Quality Scores* |
|--------------------|--|
| N (no call)        | N (no call)                                      |
| 2–9                | 6  |
| 10–19              | 15   |
| 20–24              | 22   |
| 25–29              | 27   |
| 30–34              | 33   |
| 35–39              | 37   |
| ≥ 40               | 40   |

By replacing the quality scores between 19 and 25 with a new score of 22, data storage space is conserved.

<sup>\*</sup>The mapped quality score of each bin (except "N") is subject to change depending on individual Q-tables.

#### Introduction to the BED format

- When working with BAM files, it is very common to want to examine a focused subset of the reference genome, e.g. the exons of a gene
- These subsets are commonly specified in 'BED' files (Browser Extensible Data)
   <a href="https://genome.ucsc.edu/FAQ/FAQformat.html#format1">https://genome.ucsc.edu/FAQ/FAQformat.html#format1</a>
- Many BAM manipulation tools accept regions of interest in BED format
- Basic BED format (tab separated):
  - Chromosome name, start position, end position (BED3)
  - Coordinates in BED format are 0 based:

Start position (the 0-based start coordinate)

End position (the 1-based end coordinate)

#### Introduction to the BED format

- There are several flavors of BED format: BED3, BED4, BED6, BED8, etc.
- First 3 fields always required: chr, start, stop
- Followed by up to 9 additional optional fields: name, score, strand, thickStart, thickEnd, itemRGB, blockCount, blockSizes, blockStarts

| chr7 | 127471196 | 127472363 | Pos1 | 0 | + |
|------|-----------|-----------|------|---|---|
| chr7 | 127472363 | 127473530 | Pos2 | 0 | + |
| chr7 | 127473530 | 127474697 | Pos3 | 0 | + |
| chr7 | 127474697 | 127475864 | Pos4 | 0 | + |
| chr7 | 127475864 | 127477031 | Neg1 | 0 | - |
| chr7 | 127477031 | 127478198 | Neg2 | 0 | - |
| chr7 | 127478198 | 127479365 | Neg3 | 0 | - |
| chr7 | 127479365 | 127480532 | Pos5 | 0 | + |
| chr7 | 127480532 | 127481699 | Neg4 | 0 | - |

#### Manipulation of SAM/BAM and BED files

 Several tools are used ubiquitously in sequence analysis to manipulate these files

- SAM/BAM files
  - samtools
  - bamtools
  - Picard
- BED files
  - bedtools
  - bedops



#### Common sources of confusion

Genomic coordinate systems

Genome builds

Variant representation

#### Genomic coordinates - 1 vs 0 based

| chr1    |   | Т      | Α |   | С | G |   | Т | С |   | Α |   |
|---------|---|--------|---|---|---|---|---|---|---|---|---|---|
| 1-based |   | 1<br>1 | 2 |   | 3 | 4 |   | 5 | 6 |   | 7 |   |
| 0-based | 0 |        | 1 | 2 |   | 3 | 4 | Ţ | 5 | 6 |   | 7 |

|                                      | 1-based      | 0-based      |
|--------------------------------------|--------------|--------------|
| Indicate a single nucleotide         | chr1:4-4 G   | chr1:3-4 G   |
| Indicate a range of nucleotides      | chr1:2-4 ACG | chr1:1-4 ACG |
| Indicate a single nucleotide variant | chr1:5-5 T/A | chr1:4-5 T/A |

- 1-based: Single nucleotides, variant positions, or ranges are specified directly by their corresponding nucleotide numbers
  - GFF, SAM, VCF, Ensembl browser, ...
- O-based: Single nucleotides, variant positions, or ranges are specified by the coordinates that flank them
  - BED, BAM, UCSC browser, ...

#### Reference Genome builds

#### Lift-over

Current human: GRCh38, hg38, b38


alternates: GRCh38v2\_ccdg,

GRCh38\_full\_analysis\_set\_plus\_decoy\_hla

Previous human: GRCh37, hg19, b37


**Current mouse:** mm11

Still used mouse:GRCm38, mm10



For a detailed discussion of various human reference genome flavors refer here: <a href="https://pmbio.org/module-02-inputs/0002/02/01/Reference Genome/">https://pmbio.org/module-02-inputs/0002/02/01/Reference Genome/</a>

#### Variant shifting (alignment) and parsimony/trimming



**Parsimony:** representing variant in as few nucleotides as possible without reducing the length of any allele to 0

Left (right) aligning =
shifting the start position of
a variant as far to the left
(right) as possible

### How should I sort my SAM/BAM file?

- Generally BAM files are sorted by <u>position</u>
  - This is for performance reasons
    - When sorted and indexed, arbitrary positions in a massive BAM file can be accessed rapidly
- Certain tools require a BAM sorted by <u>read name</u>
  - Usually this is when we need to easily identify both reads of a pair

45

- The insert size between two reads may be large
- In fusion detection we are interested in read pairs that map to different chromosomes

#### **Unaligned**

FASTQ (seq + quality) FASTA(seq reads only)

#### Aligned

SAM (human readable, everything, HUGE)

BAM (binary, everything)

CRAM (smaller, loss of info)

BED (smallest, text, coordinates only

size

# We are on a Coffee Break & Networking Session